

# Wide orifice pipette tips.

**NEW!**

## PROTECT DELICATE SAMPLES

Protect your delicate samples from shear force and reduce flow resistance with Starlab's wide orifice tips. These high quality pipette tips made from 100% virgin polypropylene are ideal for use in a range of applications including next generation sequencing, cell culture, and plant sciences.

- Eliminate cell fragmentation caused by mechanical shearing
- Maintain the integrity of high molecular weight DNA like gDNA or rDNA
- Ideal for pipetting viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup etc

## EASY-OPEN RACKS

No fiddly latch to catch gloves, simply press and lift the lid



[starlab.click/wideorificetips](https://starlab.click/wideorificetips)

Orifice  
inside  
diameter  
1.5 mm




200 µl

Orifice  
inside  
diameter  
1.0 mm



200 µl XL

Orifice  
inside  
diameter  
1.63 mm



1000 µl XL

## PRODUCT FEATURES

- Compatible with a wide range of pipettes
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free. Sterile product is certified endotoxin-free.
- Loose tips supplied in resealable bags
- Racked tips supplied in convenient, press-to-open, autoclavable racks

# Wide orifice pipette tips.

**NEW!**

## ORDER INFORMATION

Tips shown approx. actual size.



**200 µl** Graduated Tip / Wide Orifice

Pack Type	Sterile	Pack Size	Article No.
Bag	No	1000	<b>E1011-8000</b>
96-Tip Rack	No	10 x 96	<b>E1011-8400</b>
96-Tip Rack	Yes	10 x 96	<b>E1011-8410</b>



**200 µl** Graduated Filter Tip / Wide Orifice

96-Tip Rack	Yes	10 x 96	<b>E1011-8420</b>
-------------	-----	---------	-------------------



**200 µl XL** Graduated Tip / Wide Orifice

96-Tip Rack	No	2 x 96	<b>E1011-7400</b>
96-Tip Rack	Yes	2 x 96	<b>E1011-7410</b>



**200 µl XL** Graduated Filter Tip / Wide Orifice

96-Tip Rack	Yes	2 x 96	<b>E1011-7420</b>
-------------	-----	--------	-------------------



**1000 µl XL** Graduated Tip / Wide Orifice

Bag	No	1000	<b>E1011-9000</b>
96-Tip Rack	No	10 x 96	<b>E1011-9400</b>
96-Tip Rack	Yes	10 x 96	<b>E1011-9410</b>



**1000 µl XL** Graduated Filter Tip / Wide Orifice

96-Tip Rack	Yes	10 x 96	<b>E1011-9420</b>
-------------	-----	---------	-------------------

**Tip length.** 200 µl: 50 mm long. 200 µl XL: 91.1 mm long. 1000 µl XL: 86.4 mm long. NOTE: Filter tips are not autoclavable.

[starlab.click/wideorificetips](https://starlab.click/wideorificetips)

## PIPETTE COMPATIBILITY

ErgoOne®	HTL®
Gilson®	Nichiryo®
Biohit®	Rainin®
Brand®	Sartorius®
Capp®	Socorex®
Eppendorf®	StarPet®
Finnpipette®	StarPet® E
Hamilton®	VWR®
Hirschmann®	



Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free. Sterile product is certified endotoxin-free.