



## Thermoblock for Micro/Deepwell Plates, with Lid



#### **Product Highlights**

- Thermoblock for Starlab Mixer HC
- Automatic thermal block recognition simplifies block exchange
- For micro/deepwelll plates with SBS footprint

### **Product Applications**

Sample Preparation

Purification of plasmids, DNA, RNA

cDNA synthesis

Enzyme reactions (restriction digest, ligations, proteinase-K-digestion)

Denaturation of DNA, RNA, proteins

Labeling of DNA, RNA, proteins

Cultivation of bacteria in microcentrifuge tubes, centrifuge tubes and deepwell plates

Transformation of bacteria

Resuspension of pellets

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







### **General Data**

Art. No.	S8012-0018
Pack Size	1 Piece (1 Box × 1 Piece)
Operating temperature	4 to 35 °C
Sample capacity	1 × micro test plate / 1 × deepwell plate
Speed Range	300 – 1,400 rpm; 800 rpm when loadings > 200 g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD

5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk







# More informations about Thermoblock for Micro/Deepwell Plates, with Lid

Exchangeable thermoblock for use on the Starlab Mixer HC. Suitable for all SBS micro and deepwell plates.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD 5 Tanners Drive MK14 5BU Milton Keynes Email: infoline@starlab.co.uk

