



## Cardboard Storage Boxes 145 x 145 mm for Centrifuge Tubes



### Product Highlights

- Moisture-resistant coating
- 16-place box for 50ml tubes
- 36-place box for 15ml tubes
- Removable dividers

## Product Applications

Sample storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Dimensions (W × D × H)	145 mm × 145 mm × 80 mm
Max. tube length	140 mm
Autoclavable	No
Material	Coated Cardboard
Color	Manilla

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Cardboard Storage Boxes 145 x 145 mm for Centrifuge Tubes





Manilla cardboard boxes with lift-off lid for the storage of 15 ml or 50 ml Centrifuge Tubes. Moisture-resistant coating. Removable dividers. Max. tube height: 140 mm.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">16-Place Cardboard Storage Boxes for 50 ml Centrifuge Tubes</a> <b>Places: 16</b>  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	I9023-1650
 <a href="#">36-Place Cardboard Storage Boxes for 15 ml Centrifuge Tubes</a> <b>Places: 36</b>  <a href="#">Create PDF Data Sheet</a>	1 Piece (1 Box × 1 Piece)	I9023-1615

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)