



## Screw Cap Microtubes (Tubes only)



### Product Highlights

- Available with or without graduations
- Graduated versions marked at every 100  $\mu$ l to provide easy sample estimation
- Offered with either skirted or conical bases
- Routinely tested to 20,000 x g
- All tubes and caps are autoclavable
- Certified RNase, DNase and DNA free
- 100% medical grade polypropylene

## Product Applications

Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
RNase free	Yes
Color	Amber
g-Force max.	20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)













## More informations about Screw Cap Microtubes (Tubes only)

### Mix and match – tubes and caps!

Polypropylene screw cap tubes made from 100% medical grade polypropylene. Suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C). Available with or without graduations, the graduated versions are marked at every 100 µl to provide easy sample estimation. Offered with either skirted or conical bases. Routinely tested to 20,000 x g. Caps are available separately; all tubes and caps are autoclavable. Certified RNase, DNase and DNA free.







## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.5 ml Screw Cap Tube, Ribbed, Skirted, Amber</a> <b>Volume:</b> 0.5 ml <b>Tube Style:</b> Ribbed <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1405-2127
 <a href="#">1.5 ml Screw Cap Tube, Ribbed, Skirted, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Ribbed <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1415-2227
 <a href="#">1.5 ml Screw Cap Tube, Graduated, Conical, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Conical   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1415-2237
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Conical, Amber</a> <b>Volume:</b> 2.0 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Conical   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1420-2337
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Skirted, Amber</a> <b>Volume:</b> 2.0 ml <b>Tube Style:</b> Graduated <b>Tube Base:</b> Skirted   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Bag × 500 Tubes)	E1420-2347

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <p><a href="#">1.5 ml Screw Cap Tube, Plain, Skirted, Amber</a> <b>Volume:</b> 1.5 ml <b>Tube Style:</b> Plain <b>Tube Base:</b> Skirted</p> <p> <a href="#">Create PDF Data Sheet</a></p>	500 Tubes (1 Bag × 500 Tubes)	E1415-2247
 <p><a href="#">0.5 ml Screw Cap Tube, Plain, Skirted, Amber</a> <b>Volume:</b> 0.5 ml <b>Tube Style:</b> Plain <b>Tube Base:</b> Skirted</p> <p> <a href="#">Create PDF Data Sheet</a></p>	500 Tubes (1 Bag × 500 Tubes)	E1405-2147

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)