



## 20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, (sterile)



### Product Highlights

- Graduated at 10 µl
- UltraPoint® thin-wall orifice helps prevent sample hang up
- Highly effective HDPE filter that contains no additive to interfere with your samples
- Universal fit. One tip fits all! The tip collar is engineered to fit all popular pipettes
- Certified RNase, DNase, DNA and Pyrogen free

## Product Applications

DNA applications (e.g. PCR)

RNA applications (e.g. Gene expression analysis)

Protein applications (e.g. Antibody Research)

Cell Culture applications (e.g. Media)

Applications with radio actives

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)



## General Data

Art. No.	See variations
Pack Size	96 Tips (96 Tips)
Sterile	Yes
Autoclavable	No
Filter tip	Filter
Volume range	0.1 - 20 µl
Volume	20 µl
Tip Type	UltraPoint®
Material	Polypropylene
Graduations	Yes
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about 20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, (sterile)

### Nothing goes through!

Protect your samples against unwanted contamination by using TipOne® filter tips. The TipOne® filter is one of the best on the market and protects against aerosol cross-contamination without trapping your valuable samples. The filters contain no additives to interfere with samples.








During pipetting it is not possible to see with the eye if fine aerosols are going into your pipette. If they are, those aerosols can contaminate other samples at a later date.

TipOne® filter tips are proven to block >99% of aerosols and are among the best filters on the market. The results, shown right, speak for themselves!

*\*The comparative tests were conducted in an independent test institute with filter tips from other leading manufacturers.*



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Refills (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Refills × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Refills × 96 Tips)	S1123-1710
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Refills (sterile), Case</a> <b>Pack Size:</b> 7,680 Tips (80 Refills × 96 Tips)	7,680 Tips (80 Refills × 96 Tips)	S1123-1710-C
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile)</a> <b>Pack Size:</b> 960 Tips (10 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	960 Tips (10 Racks × 96 Tips)	S1123-1810
 <a href="#">20 µl UltraPoint® graduated TipOne® Filter Tip, Natural, Racks (sterile), Case</a> <b>Pack Size:</b> 7,680 Tips (80 Racks × 96 Tips)  <a href="#">Create PDF Data Sheet</a>	7,680 Tips (80 Racks × 96 Tips)	S1123-1810-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)