



Media Bottle, PS, 1000 ml (sterile)



Product Highlights

- 1000 ml PS Media Bottle
- Compatible with all common Bottle Top Filters (GL 45 compatible)
- Wide base for stability
- Temperature range: -20 °C to +50 °C
- Sterile SAL10⁻⁶ (ISO1137) / 3-year shelf life
- Certified RNase, DNase free and non-pyrogenic
- Wide bottle top for fast and efficient use
- Ergonomically designed for easier handling to help prevent spills.

Product Applications

Sterile filtration
Media preparation
Media storage

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



General Data

Art. No.	CC8221-6000
Pack Size	24 Pcs. (1 Box × 24 Pcs.)
Sterile	Yes
Autoclavable	No
Volume	1000 ml
Material	Polystyrene
DNase free	Yes
Free of endotoxins	Yes
Free of heavy metals	Yes
Pyrogen free	Yes
RNase free	Yes
Non-cytotoxic	Yes
System	Media Bottle

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



More informations about Media Bottle, PS, 1000 ml (sterile)

OPEN WIDE!






Media Bottles with a wide base for stability and a wide, easily accessible bottle top for efficient use. The reservoir bottle caps have large knurls for easy loosening or tightening. Each product is manufactured in a Class 100,000 cleanroom environment and is printed with the batch number and expiry date for easy traceability.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk



All Variations

PRODUCT NAME		PACKAGING SIZE	ART. NO.
	Media Bottle, PS, 150 ml (sterile) Volume: 150 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8221-3000
	 Create PDF Data Sheet		
	Media Bottle, PS, 250 ml (sterile) Volume: 250 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8221-4000
	 Create PDF Data Sheet		
	Media Bottle, PS, 500 ml (sterile) Volume: 500 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8221-5000
	 Create PDF Data Sheet		
	Media Bottle, PS, 1000 ml (sterile) Volume: 1000 ml	24 Pcs. (1 Box × 24 Pcs.)	CC8221-6000
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab (UK) LTD
5 Tanners Drive
MK14 5BU Milton Keynes
Email: info@starlab.co.uk