



## Screw Cap Microtubes



### Product Highlights

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## General Data

Art. No.	See variations
Autoclavable	Yes
Material	Polypropylene
DNA free	Yes
DNase free	Yes
Free of endotoxins	Yes
PCR inhibitor free	Yes
RNase free	Yes
Color	Natural
g-Force max.	20,000 x g

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



## More informations about Screw Cap Microtubes

### Mix and match – tubes and caps!











Polypropylene screw cap tubes that are suitable for many applications from sample storage in the gas phase of liquid nitrogen to freezer, through boiling to autoclave (121 °C).

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [info@starlab.co.uk](mailto:info@starlab.co.uk)













## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">0.5 ml Screw Cap Tube, Ribbed, Skirted, Natural</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Ribbed   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1405-2120
 <a href="#">0.5 ml Screw Cap Tube, Plain, Conical, Natural</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Plain   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1405-2130
 <a href="#">0.5 ml Screw Cap Tube, Plain, Skirted, Natural</a> <b>Volume:</b> 0.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Plain   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1405-2140
 <a href="#">1.5 ml Screw Cap Tube, Ribbed, Conical, Natural</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Ribbed   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1415-2210
 <a href="#">1.5 ml Screw Cap Tube, Ribbed, Skirted, Natural</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Ribbed   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1415-2220

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">1.5 ml Screw Cap Tube, Graduated, Conical, Natural</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1415-2230
 <a href="#">1.5 ml Screw Cap Tube, Plain, Skirted, Natural</a> <b>Volume:</b> 1.5 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Plain   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1415-2240
 <a href="#">2.0 ml Screw Cap Tube, Ribbed, Skirted, Natural</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Ribbed   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1420-2320
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Conical, Natural</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Conical <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1420-2330
 <a href="#">2.0 ml Screw Cap Tube, Graduated, Skirted, Natural</a> <b>Volume:</b> 2.0 ml <b>Tube Base:</b> Skirted <b>Tube Style:</b> Graduated   <a href="#">Create PDF Data Sheet</a>	500 Tubes (1 Pack × 500 Tubes)	E1420-2340

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab (UK) LTD**  
5 Tanners Drive  
MK14 5BU Milton Keynes  
Email: [infoline@starlab.co.uk](mailto:infoline@starlab.co.uk)