



200 µl TipOne® Tip, Yellow, Refills (non-sterile)



Product Highlights

- Universal tip fit
- High quality virgin polypropylene tips, with extremely smooth surfaces for reduced sample retention
- Certified RNase, DNase, DNA and Pyrogen free
- Refills dramatically reduce plastic waste compared to pre-racked tips
- Refill top and bottom shells made from 100% recycled plastic
- All components, excluding tips and shrink wrap, are recyclable

Product Applications

Pipetting liquids

Dispensing liquids

Mixing liquids

Filling plates and reaction vessels

Phase extraction and removal of supernatants

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	S1111-0706
Pack Size	960 Tips (10 Refills × 96 Tips)
Color Code Wafer	Yellow
Sterile	No
Autoclavable	Yes
Filter tip	Non-filter
Volume range	1.0 - 200 µl
Volume	200 µl
Tip Type	standard
Material	Polypropylene
Graduations	No
Color	Yellow
Sustainability attributes	Refills dramatically reduce plastic waste compared to pre-racked tips. Refill top and bottom shells made from 100% recycled plastic. All components, excluding tips and film wrap, are recyclable

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



For use with

Abimed: BioPette-E, Discovery Comfort, LABMATE Optima

Biohit: eLine, ePet, mLine, Picus, Proline

BRAND: Transferpette, Transferpette Electronic, Transferpette Pro, Transferpette S

Capp: CappAero Comfort, CappAero Multi, CappAero Universal, Capptronic Multi, Capptronic Single, ecopipette, Singlechannel

Eppendorf: Reference, Reference 2, Research, Research plus, Research pro, Xplorer

Gilson: Pipetman Concept, Pipetman L, Pipetman Neo, Pipetman P, Pipetman Ultra

Hamilton: SoftGrip Pipette, SofTouch

HTL: Clinipet+, Discovery Comfort, Discovery Pro, Labmate Pro, Optipette

Nichiryo: 00-NP7, 50F, 50V, NLE, NPLO, NPX

Rainin: Classic, EDP1, EDP-3 Plus, Pipet-Lite, Pipet-Plus

Sartorius: eLine, ePet, mLine, Picus, Picus NxT, Proline, Proline Plus, Tacta

Socorex: Acura 815 manual, Acura 825, Acura 855, Acura electro 925, Acura electro 955, Calibra digital

STARLAB: ErgoOne, ErgoOne E, MicroOne, StarPet, StarPet-E

Thermo Fisher Scientific Inc.: Finn Focus, Finn Focus Micro, Finn Focus Univ., Finn pipette Digital, Finn pipette F1/F2, Finn pipette Novus, Finn pipette Novus Turquoise, Finn pipette Novus Yellow

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about 200 µl TipOne® Tip, Yellow, Refills (non-sterile)

High quality, reliable pipette tips with a universal tip fit – TipOne® has it all! If you are looking for a sustainable pipette tip system with accuracy and precision you can trust, TipOne® is the tip for you.

TipOne® refills are a quick and easy way to obtain economical racked tips, and they save on plastic use!

Reduce single-use plastic!













Compared to pre-racked tips, using TipOne® refills dramatically reduces the amount of single-use plastic going through your laboratory. Quick and easy to load, unwrap and click in place into an empty TipOne® rack. The refills are supplied in two separate towers of 5x96 tips, so you can rack up 5 or 10 at a time. Tip information is printed on every wafer, providing full traceability when racked. Once racked, and if needed, the refills can be autoclaved.

Join our journey to sustainable pipetting

The clear, polypropylene system components are made from 100% recycled plastic. TipOne® tip racks are built to last and can be reused again and again with these cost-effective refills. All components of TipOne® refills, excluding tips and shrink wrap, are recyclable. Sign up for our TipOne® recycling service!











All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 200 µl TipOne® Tip, Yellow, Racks (non-sterile) Sterile: No Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-0806
 200 µl TipOne® Tip, Yellow, Racks (non-sterile), Case Sterile: No Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-0806-C
 200 µl TipOne® Tip, Yellow, Racks (sterile) Sterile: Yes Pack Size: 960 Tips (10 Racks × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Racks × 96 Tips)	S1111-0816
 200 µl TipOne® Tip, Yellow, Racks (sterile), Case Sterile: Yes Pack Size: 7,680 Tips (80 Racks × 96 Tips)  Create PDF Data Sheet	7,680 Tips (80 Racks × 96 Tips)	S1111-0816-C
 200 µl TipOne® Tip, Yellow, Refills (non-sterile) Sterile: No Pack Size: 960 Tips (10 Refills × 96 Tips)  Create PDF Data Sheet	960 Tips (10 Refills × 96 Tips)	S1111-0706
 200 µl TipOne® Tip, Yellow, Refills (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (100 Refills × 96 Tips)  Create PDF Data Sheet	9,600 Tips (100 Refills × 96 Tips)	S1111-0706-C

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME		PACKAGING SIZE	ART. NO.
	200 µl TipOne® Tip, Yellow, Stack rack (non-sterile) Sterile: No Pack Size: 960 Tips (1 Stack rack × 960 Tips)	960 Tips (1 Stack rack × 960 Tips)	S1111-0206
	 Create PDF Data Sheet		
	200 µl TipOne® Tip, Yellow, Stack racks (non-sterile), Case Sterile: No Pack Size: 9,600 Tips (10 Stack racks × 960 Tips)	9,600 Tips (10 Stack racks × 960 Tips)	S1111-0206-C
	 Create PDF Data Sheet		
	200 µl TipOne® Tip, Yellow, Bags (non-sterile) Sterile: No Pack Size: 1,000 Tips (1 Bags × 1000 Tips)	1,000 Tips (1 Bags × 1000 Tips)	S1111-0006
	 Create PDF Data Sheet		
	200 µl TipOne® Tip, Yellow, Bags (non-sterile), Case Sterile: No Pack Size: 25,000 Tips (25 Bags × 1000 Tips)	25,000 Tips (25 Bags × 1000 Tips)	S1111-0006-C
	 Create PDF Data Sheet		

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltingbaum 38
22143 Hamburg
Email: info@starlab.de