




Large Volume Tips 5 ml



Product Highlights

- From 1200 µl tips to 10 ml tips
- Made from polypropylene
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



General Data

Art. No.	See variations
Autoclavable	Yes
Filter tip	Non-filter
Volume	5 ml
Material	Polypropylene

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



More informations about Large Volume Tips 5 ml













Large Volume Tips

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de















All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 5 ml Biohit-Style Tip (non-sterile) Design Tips: Biohit-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet	250 Tips (1 Bag × 250 Tips)	E1009-7000
 5 ml Biohit-Style-XL Tip (non-sterile) Design Tips: Biohit-Style-XL Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet	250 Tips (1 Bag × 250 Tips)	E1009-7001
 5,000 µl Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet	250 Tips (1 Bag × 250 Tips)	I1050-0000
 5 ml Gilson-Style Graduated Macro Tip (non-sterile) Design Tips: Gilson-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1050-0700
 5 ml Gilson-Style Graduated Macro Tip (sterile) Design Tips: Gilson-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1050-0710
 5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet	250 Tips (1 Bag × 250 Tips)	I1053-0000

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de



PRODUCT NAME	PACKAGING SIZE	ART. NO.
 5 ml Eppendorf-Style Graduated Macro Tip (non-sterile) Design Tips: Eppendorf-Style Sterile: No Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1053-0700
 5 ml Eppendorf-Style Graduated Macro Tip (sterile) Design Tips: Eppendorf-Style Sterile: Yes Pack Size: 500 Tips (10 Racks × 50 Tips)  Create PDF Data Sheet	500 Tips (10 Racks × 50 Tips)	I1053-0710
 5 ml Finn-Style Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 250 Tips (1 Bag × 250 Tips)  Create PDF Data Sheet	250 Tips (1 Bag × 250 Tips)	I1009-5000
 5 ml Finn-Style Pipette Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 200 Tips (1 Bag × 200 Tips)  Create PDF Data Sheet	200 Tips (1 Bag × 200 Tips)	E1070-2596
 5 ml Finn-Style Pipette Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 1,000 Tips (5 Bags × 200 Tips)  Create PDF Data Sheet	1,000 Tips (5 Bags × 200 Tips)	E1070-2601
 5 ml Finn-Style Pipette Tip (non-sterile) Design Tips: Finn-Style Sterile: No Pack Size: 28 Tips (1 Rack × 28 Tips)  Create PDF Data Sheet	28 Tips (1 Rack × 28 Tips)	E1070-2607

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

Starlab International GmbH
Neuer Höltigbaum 38
22143 Hamburg
Email: info@starlab.de