




## 10 µl Round Gel-Loading Tips



### Product Highlights

- Flexible, long and narrow
- Loading samples onto gels (eg. polyacrylamide or agarose)
- Handling of small volumes
- Compatible with a wide range of pipettes
- Different pack types available
- Certified

## Product Applications

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## General Data

Art. No.	See variations
Filter tip	Non-filter
Volume	10 µl
Tip Type	round gel-loading
Color	Natural

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## More informations about 10 µl Round Gel-Loading Tips











All tips are certified RNase and DNase free. Packaging of these tips differs from the TipOne® range.

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)



## All Variations

PRODUCT NAME	PACKAGING SIZE	ART. NO.
 <a href="#">10 µl Round Gel-Loading Tip, Loose (non-sterile)</a> <b>Sterile:</b> No <b>Length:</b> 56.62 mm   <a href="#">Create PDF Data Sheet</a>	1,000 Tips (1 Bag × 1000 Tips)	I1010-3000
 <a href="#">10 µl Round Gel-Loading Tip, Racked (Sterile)</a> <b>Sterile:</b> Yes <b>Length:</b> 56.62 mm   <a href="#">Create PDF Data Sheet</a>	1,020 Tips (5 Racks × 204 Tips)	I1010-3610
 <a href="#">10 µl Round Gel-Loading Tip, Racked (non-sterile)</a> <b>Sterile:</b> No <b>Length:</b> 71.25 mm   <a href="#">Create PDF Data Sheet</a>	1,632 Tips (8 Racks × 204 Tips)	I1012-3600
 <a href="#">10 µl Round Gel-Loading Tip, Racked (Sterile)</a> <b>Sterile:</b> Yes <b>Length:</b> 71.25 mm   <a href="#">Create PDF Data Sheet</a>	1,632 Tips (8 Racks × 204 Tips)	I1012-3610
 <a href="#">10 µl Round Gel-Loading Tip, Racked (non-sterile)</a> <b>Sterile:</b> No <b>Length:</b> 68 mm   <a href="#">Create PDF Data Sheet</a>	1,020 Tips (5 Racks × 204 Tips)	I1022-8950

STARLAB reserves the right to make changes at any time and without prior notice. The content and design of this PDF are protected by national and international copyright law and are the property of STARLAB International GmbH. Any duplication, editing, distribution and any kind of use and utilization of this PDF content in electronic systems, online media and / or libraries or similar databases requires the prior consent of STARLAB International GmbH.

**Starlab International GmbH**  
Neuer Höltigbaum 38  
22143 Hamburg  
Email: [info@starlab.de](mailto:info@starlab.de)